



Surgical perfection. For life.

STILLE Medstone URO/ GYN Urology/ Cysto/ Brachytherapy/ Gynecology

High Quality Hospital Grade C-Arm Tables for Urology and Gynecology

The STILLE Medstone c-arm tables contribute to enhancing efficiency while representing an affordable solution with quality and reliability provided by Swedish technology.

STILLE Medstone URO Surgical Imaging Tables

Max. 4 motorized movements options available:

- Height adjustment 25.98"-39.76" (66-101cm)
- Trendelenburg $\pm 15^\circ$
- Longitudinal travel 15.75" (40cm)
- Lateral travel 6" (15cm) (only Medstone4 URO)



Photo shows Medstone4 URO (500-1028) with leg holders 505-1126

STILLE Medstone4 URO Standard Configuration and Specifications

- Carbon fiber table top 24"x48" (61x122cm) with 1mm Aluminium Equivalence
- Imaging area 24"x24" (61cmx61cm)
- Weight capacity of the table top 500lb (227kg)
- 230Vac, 50 Hz, 3.4A (120Vac, 60Hz, 6.5A) with 230Vac (120Vac) style line cord
- Four hospital grade locking casters
- Infectious control covers
- 2" (5cm) table mattress
- Wired hand control
- Wired foot control with pedestal to place the foot control on the table base
- One flexible uro-catcher drain bag mounting hoop
- One starter package of 20 disposable drain bags

Model	4 URO
Part number	500-1028
Table top	24"x48" (61cm x 122cm)
Height adjustment	25.98"-39.76" (66-101cm)
Trendelenburg	$\pm 15^\circ$
Lateral tilt	x
Longitudinal travel	40cm
Lateral travel	15cm
Auto horizontal and center home positioning	✓
Type of cover	Infectious Control

Specification tolerance: $\pm 1^\circ$, $\pm .5"$ (1.27cm)
Medstone URO tables are also available with 2 or 3 movements, please contact us for further information.



Optional pair of gas spring assisted leg holders 505-1114

Optional Accessories

Part number	Description
505-1126	Leg holders Goepel type/ pair with attachment clamps
505-1114	Leg holders/ pair with gas spring support and quick lock release handle
505-1125	Radiolucent leg extension plate 24"x28" (61cmx71cm), 16.7lbs (7.6kg) incl. pad
505-1170	Non-radiolucent leg extension plate 24"x25" (61x64cm), 9.5lbs (4.3kg) incl. pad
514-50-3SR025	Universal adjustable arm board with pad and clamp
505-1127	10 non-sterile drain bags
505-1120	Replacement uro-catcher drain bag mounting loop

STILLE – Surgical perfection since 1841

STILLE is one of the oldest medical device companies in the world founded at Karolinska University Hospital in Sweden back in 1841. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas – Surgical Instruments and Surgical imaging tables.

Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery. Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2 surgical imaging table has the highest translucency in the world, and is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures, enhancing crystal imaging while reducing radiation exposure as well as increasing efficiency in the OR.

The STILLE Medstone line of surgical imaging tables range from floating to fully motorized tables for vascular, electrophysiology, gastro (ERCP), pain management and urology to suit your needs, with reliability and optimal performance at an affordable cost at the highest standards of Swedish technology.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
sales@stille.se

www.stille.se

STILLE Surgical Inc.

915 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614
sales.us@stille.se