

Optional Accessories

Features:

All table accessories are made of SUS304 stainless steel. All the surfaces are easy to clean and disinfect. With 6 standard and optional accessories, the surgical bed can support a multitude of body positions and procedures, meeting the requirements of most surgeries.

 <p>121000347/121000348/ Foot plate (notched clamp)</p>	 <p>121000349/ High position arm holder</p>	 <p>121000351/ Connection bracket</p>
 <p>* 121000350/121000351/ Head rest(connection bracket)</p>	 <p>121000352/121000351/ Head plate(connection bracket)</p>	 <p>121000353/ Hand operating table</p>
 <p>121000354/ Instrument tray</p>	 <p>121000355/ Iliac support</p>	 <p>121000356/ Intervertebral pad</p>
 <p>121000357/ Drainage bowl</p>	 <p>* 121000358/ Orthopedics traction</p>	 <p>121000360/ Infusion frame</p>

Orthopedics Features

The orthopedics traction is designed based on EU's requirements, meeting the requirements of most orthopedics surgeries. All parts are made of SUS304 stainless steel with high strength. Therefore, the surgical table is very stable during surgeries and is X-ray and C-arm translucent.

Standard Accessories

 <p>121000248/ Round-hole clamp</p>	 <p>121000247/ Anesthesia screen</p>	 <p>121000251/ Body strap</p>
--	---	--



Gel Pad

Functions & Features:

The gel pad which made from viscoelastic polymer gel, offers stability and comfort during positioning. It has uniform dispersion of pressure which can avoid bedsores. It is easy to clean and X-ray translucent. It can support surgery with working environment temperature from -10°C to 50°C. It is also biocompatible with human skin and tissues. All these features make the gel pads cost-effective and durable.

