



## E-Z Lift Beach Chair

*Position patients for shoulder surgery with the lightweight shoulder arthroscopy beach chair from SchureMed*

Position patients for shoulder surgery with this lightweight and cost-effective surgical beach chair that adjusts from 0-90° with a powerful piston system that assists O.R. staff when lifting the patient into position. A safe and secure option to our fully-powered shoulder surgery chair system.

### Features

- Significantly reduces effort by O.R. staff when lifting patients over other manual systems
- Double-ball joint head system safely secures patient's head
- 500 lb. surgical patient weight capacity
- Articulation range 0-90°
- Sliding head rest conforms to every patient
- Five-piece Deluxe Pad Set & Lateral Brace included
- 28"W x 30"D x 47"H, 32 lbs



#800-0142 E-Z Lift Beach Chair,  
#800-0142-EU E-Z Lift Beach Chair, EU

### Required Accessories: Select One Model

- #800-0164 Disp. Face Masks, 12/cs
- #800-0039 Disp. Head & Chin Straps, 12/cs
- #508-1345 Reusable Head & Chin Straps, set

### Optional Accessories

- #800-0074 Shoulder Chair Dolly
- #800-0050 Multi-Axis Arm Positioner
- #800-0134 Schure Socket XL

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)  
[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351