



## Uro-Trapper System


*Lock fluid drainage device into place w/easy-mount clamps*

Patient positioning is never compromised nor is the surgeon's access to the patient with this flexible, surgical fluid containment system that attaches to any general surgical table.

### **Features**

- *Contains fluids during hysteroscopy, cystoscopy, urologic and oncologic surgical procedures*
- *Flexible drape support provides unobstructed surgeon access to surgical site*
- *Once compressed by surgeon's weight during procedure, the support springs back to allow optimal fluid drainage*
- *Mounts effortlessly to any surgical O.R. table with integral clamps*
- *Clamps are included*

**#800-0096** Uro-Trapper System



**Phone** 888.724.8763 • 781.982.7000 | **Fax** 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351