



Great White® Platinum Stirrups

Gynecology, laparoscopy & urology procedures are safer than ever

The Great White® Platinum Stirrups provide safe and quick positioning assistance for improved surgical site access while offering exceptionally precise placement.

Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -30° low to 75° high; Adduction -9° and abduction +25°
- The raised lateral fin boots effectively minimize the possibility of nerve damage
- Included Platinum Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-PL Great White® Platinum Stirrups

Required Accessory: Purchase Two SpringLoc Clamp

- | | |
|-------------------|-----------------------|
| #800-0338 US | #800-0338-EU EU |
| #800-0338-JPN JPN | #800-0338-SWISS SWISS |

Optional Accessories

- #508-1415 Replacement Stirrup Boots Pads w/Fins
- #800-0074-S Great White® Stirrup Dolly
- #800-0284 Stirrup Wall Rack
- #800-0160 Clear View Leg Drapes

Improved release mechanism in handles

New Improved Design!

600 lb. (272kg)
Patient Weight Capacity



Redesigned mechanical ball joint



Increase your surgeons' O.R. performance

Patent Pending

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351